

UNICORN

Computer specifications

	UNICORN™ Client	Database Server	Workstation installation	E-License Server
Min. free disk space	16 GB	10 GB	27 GB	500 MB
Min. available RAM	4 GB	4 GB	8 GB	2 GB
Disc format	NTFS	NTFS	NTFS	NTFS
Architecture	Quad-core processor or 4 logical processor	Quad-core processor or 4 logical processor	Quad-core processor or 4 logical processor	Multi-core processor



cytiva.com/unicorn

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

UNICORN is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

UNICORN © 2020 Cytiva

Any use of UNICORN is subject to Cytiva Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

29283583 AD V:3 11/2020